

## Electrophoretic NMR Instrumentation and Methodology

Extending the capability of conventional NMR instruments

### eNMR experimental setup



Main unit with encoder, high voltage amplifier, and the safety protection system

### Experimental setup. Block diagram



The block diagram of the eNMR assembly

## Methodology



Detection of displacement yields the electrophoretic mobility. The combination of the self-diffusion coefficient and the electrophoretic mobility yields the effective charge, a measure of ion association.<sup>1</sup>

### Methodology

Electrophoretic NMR (eNMR) - based on pulsed field gradients, akin to diffusion NMR

A double-stimulated-echo pulse sequence for eNMR



0

The signal is phase modulated:

$$\frac{s}{s_o} \propto \exp(-\gamma^2 g^2 \delta^2 D\Delta) \exp(i\phi), \quad \phi = \gamma g \delta \mu E\Delta$$
$$\mu = \frac{\phi}{\gamma \delta g \Delta_E E}$$

Resulting series of spectra

### **Experiment and Operation**



The electrophoretic mobility and the effective charge can be obtained by recording the variation of spectral phase by increasing electric field

## Experimental artefacts and the way to solve it

Some experimental artifacts may cause big signal attenuation and the artificial phase shift:

- Thermal convection in highly conductive samples
- Electroosmosis
- Gas buble production due to electrolysis





### Experimental artefacts and the way to solve it



### **Experimental artefacts and solutions**



(a)



- 100 mM TMABr in  $D_2O$
- PNIPAM with 200 mM NaSCN in  $D_2O$ solution.

- a) The porous plug retains the electrode reaction products out of the sensitive region. TMABr
- b) Gas bubbles are confined by the plug at the bottom. PNIPAM sample

#### Attenuation of the signal with and without the



The symmetric construction of the eNMR cell and the filter plug helps <sup>2</sup>. The current stabilization mode in eNMR unit is required

# Application. 2D electrophoretic mobility ordered spectroscopy (2D-MOSY) <sup>3</sup>

A mixture of L-Lysine (I), L-Serine (s) and L-Aspartic acid (a)



## eNMR 2D-MOSY for analytical applications <sup>3</sup>

Analysis of components of a drug





### **Application examples**

Investigation of zwitterionic behaviour in different pH



### Specifications of eNMR unit

| Output voltage:                              | 0 to ±1000 V                        |
|--|-------------------------------------|
| Digital To Analogue converter size: 2×12 bit |                                     |
| Output current                               |                                     |
| at ±1000 V:                                  | 0 to $\pm$ 50 mA                    |
| at ±500 V:                                   | 0 to ±200 mA                        |
| Output power:                                |                                     |
| Peak power                                   | 300 W                               |
| Mean power                                   | 30 W                                |
| Minimum/maximum pulse length:                | 500 μs / 30 s                       |
| Duty cycle:                                  | 30 %                                |
| Output pulse shapes:                         | Rectangular                         |
| Slew rate:                                   | Greater then 25 V / $\mu s$         |
| Settling time (to 2% ):                      | Less than 200 $\mu s$ for 2 kV step |
| Stability                                    |                                     |
| Drift with time                              | Less than 100 ppm/hr, noncur        |
| Drift with temperature                       | Less than 300 ppm/°C                |
| RF Filters                                   | 10 MHz low pass                     |
| Power consumption:                           | 80 W for 220 V AC                   |
| Dimensions of the main unit:                 | 430×340×90 (fits in standard        |
| Weight:                                      | 7 kg                                |
|  |                                     |



han 100 ppm/hr, noncumulative han 300 ppm/°C Hz low pass for 220 V AC 340×90 (fits in standard 19" rack)

### Key features

- Well-established scientific background and proven technology
- Operates as add-on for any conventional NMR spectrometer and probe; requires no additional hardware or software
- Advanced sample cell and RF filter system
- Straightforward embedding of high voltage pulses in conventional NMR pulse programs
- Probe and user protection system

### Selected applications of electrophoretic NMR

**Physical chemistry** – ion pairing and association in simple and complex (polyelectrolytes) ionic mixtures.

Batteries and fuel cells - chemically selective measurement of ionic migration.

**Biochemistry** – biomolecular charge and association.

Analytical chemistry – electrophoretic analysis of complex ionic mixtures

**Pharmaceutical chemistry** – release and association of charged drugs

**Metallorganic chemistry** – the structure of supramolecular complexes from the observed charge

### Conclusions

eNMR together with diffusion NMR provide a powerful tool to study the ionic solutions and complexation/binding of ions with macromolecules in solution. It is especially good at detecting weak interactions

Two dimensional eNMR (MOSY) is superior than DOSY in conventional chemical analysis and has great potential to be used in drug discovery and metabolomics

### References



### Selected literature



1. Peter Stilbs. Diffusion and Electrophoretic NMR, Series: De Gruyter STEM 2019, ISBN : 3-11-055165-9 ISBN : 3-11-055153-5

- Yuan Fang, Pavel V Yushmanov, Istvan Furo, Improved accuracy and precision in electrophoretic NMR experiments. Current control and sample cell design; July 2020, Journal of Magnetic Resonance
- 3. Y. Fang, P. V. Yushmanov, and I. Furó, Assessing 2D electrophoretic mobility spectroscopy (2D MOSY) for analytical applications, *Magn. Reson. Chem.* **55** DOI: 10.1002/mrc.4558 (2017).
- 4. M. Giesecke, F. Hallberg, Y. Fang, P. Stilbs, and I. Furó, Binding of monovalent and multivalent metal cations to polyethylene oxide in methanol probed by electrophoretic and diffusion NMR, *J. Phys. Chem. B* **120** 10358–10366 (2016).
- 5. M. Bielejewski, M. Giesecke and I. Furó, On electrophoretic NMR. Exploring high conductivity samples, *J. Magn. Reson.* **243** 17-24 (2014).
- L. Patel, O. Mansour, M. Crossman, P. Griffiths, Electrophoretic NMR characterization of charged side chain cationic polyelectrolytes and their interaction with the anionic surfactant, sodium dodecyl sulfate, *Langmuir* 2019, 35, 28, 9233-9238